

## Sample protocol for obtaining assimilation-temperature (AT) response curves using the Fast Assimilation-Temperature Response (FAsTeR) method with the LI-6800.

Written by Josef Garen, last modified 5 September 2023.

1. Open a log file.
  - a. Under Log Files > Logging Options: Select “Also log to Excel file” and “Use additional averaging time” set to 1 sec. We use minimal signal averaging time because we are measuring in non-equilibrium conditions, but using 1 sec averaging time substantially reduces noise compared with no signal averaging.
  - b. Under Log Files > Matching Options: Select “Never match” for both CO<sub>2</sub> and H<sub>2</sub>O.
  - c. Under Log Files > Fluorometer Options: Select “0: Nothing” from the drop-down menu.
  - d. Open a new file to store data in. Under Log Files > Open a log file: Select a folder/make a folder, and press “New File” to make a new logfile. Use desired naming convention for the logfile.
2. Set Environment controls
  - a. Flow: Set Flow to 600  $\mu\text{mol s}^{-1}$  and  $\Delta P$  to 0.1 kPa.
  - b. H<sub>2</sub>O: Set H<sub>2</sub>O on; set RH to 35%.
  - c. CO<sub>2</sub>: Set CO<sub>2</sub> injector to on; set reference CO<sub>2</sub>R to 420ppm.
  - d. Under Environment > Fan: Set mixing fan to on; set fan speed to 10,000 rpm.
  - e. Under Environment > Temperature: Set Temperature on; Set Txchg = 0 °C
  - f. Under Environment > Light > Fluorometer: Set control mode to “Setpoint”; Set Setpoint to saturating PAR.
3. If Tair <= ambient, insert a leaf and close the chamber, ensuring that the leaf is in contact with the leaf thermocouple.
4. Under Environment > Flow, check that Leak % < 2-3%. (Do not use chamber leak test.) If Leak % is too high, fix leaks by reclamping or sealing around stems/midribs/veins with putty. You may also be able to reduce leak % by manually pressing on the chamber, but this may damage fragile leaves. Once Leak % < 2-3%, proceed.
5. Under Measurements: Monitor A, gsw, Txchg, and Tleaf for a **minimum of 20 minutes**, and then ensure stability using graph zoom functions (flat lines for 5 minutes). Some species may take more than 45 minutes to stabilize.
6. Ensure LI-COR is able to hit all setpoints: Go to the Environment tab, and check each environmental variable that you have set (e.g. RH, CO<sub>2</sub>R, Txchg) and ensure that the actual value is the same as the setpoint/requested value (or very close).
7. Check  $\Delta\text{CO}_2 < -8 \mu\text{mol s}^{-1}$ . If not, reduce flow rate until  $\Delta\text{CO}_2 < -8 \mu\text{mol s}^{-1}$  is achieved, to a minimum of  $\sim 250 \mu\text{mol s}^{-1}$  (flow in both sample and reference IRGAs must exceed  $200 \mu\text{mol s}^{-1}$ ). If  $\Delta\text{CO}_2 < -8 \mu\text{mol s}^{-1}$ , increase flow rate to bring  $\Delta\text{CO}_2$  close to (but below)  $-8 \mu\text{mol s}^{-1}$ . High flow rates are desirable in this method because they reduce the mixing and transit time for gas in the chamber, which reduces the temporal “mismatch” between the CO<sub>2</sub> measurements and the Tleaf measurement. A flow in excess of  $1000 \mu\text{mol s}^{-1}$  is likely unnecessary.
8. When Txchg is stable at its minimum temperature, set Txchg to its current value.
9. Configure temperature ramp.

- a. Under Environment > auto controls, select “Txchg”. Set F(t) to linear. Set F(0) to the current Txchg value and set F(end) to the desired maximum Txchg temp (usually 40 degrees above minimum).
  - b. Set Period based on desired ramping speed.
    - a. If flow  $\geq 600$ , use 2 degree/min, which is a period of 20 minutes for a 40 degree span. This is the fastest ramping speed that is reliably achievable by the machine.
    - b. If using a lower flow rate, use 1.5 degree/min, which is a period of 30 min for a 40 degree span. This helps reduce “mismatch” between measured temperatures and gas exchange rates.
  - c. **Press the right-facing arrow on the bottom left.** This ensures that you start the ramp from the bottom.
10. Configure autologging program
- a. Under the Auto Programs tab, select AutoLog from the drop-down box.
  - b. Select desired logging interval. Logging every 2 seconds is suggested – you can log every 1 s, but this causes the machine to bog down.
  - c. Select total time to equal or exceed the Period selected in Step 9. – If using a fast ramping speed, add an extra 5 minutes (total 25 min), and if using a slower ramping speed, add an extra 3 min (total 33 min) to the ramping period to ensure that a complete curve is logged, since often Tleaf lags behind Txchg by a few minutes.
11. Monitor Txchg, A, and gsw using charts (flat line for 5 minutes) and stability metrics (line e, “Stable” = “3/3”). When stability is achieved, proceed.
12. **Match.** Wait for stability.
13. Begin data collection. Go to Environment > Auto Controls and press **start**. Immediately go to Auto Programs and press **start**. The machine will begin ramping the temperature and logging automatically.
14. Wait while data collection is carried out.
15. **Match again at the end and log once** to see difference in match offsets.
17. Close Log file.